

Water vapor retrieval from Zugspitze FTIR measurements

Outline

- Zugspitze/Garmisch instrumentation
- Interference errors and micro window selection
- Zugspitze FTIR profile retrieval
- Strategy for validation of FTIR total columns by sondes
- Optimized FTIR total column retrieval and its validation
- Summary and Outlook

R. Sussmann, T. Borsdorff, and C. Camy-Peyret

Research Center Karlsruhe IMK-IFU, Garmisch-Partenkirchen, Ralf Sussmann

Water vapor retrieval from Zugspitze FTIR measurements

Triple NDACC Primary Station:
FTIR, Aerosol, UV.

Permanent Ground-Truthing Facility
Zugspitze/Garmisch according to the WMO
requirements.

IMK-IFU Working Group
„Variability and Trends“

Scientists

R. Sussmann
H.E. Scheel
T. Trickl
H. Vogelmann


Engineers

H. Giehl
M. Rettinger

PhD students

T. Borsdorff
F. Forster


MAPS,
CRISTA,
MOPITT,
SAGE,
GOME,
AIRS,
SCIAMACHY,
ACE,
IASI,
CALIPSO,
OCO,
TCCON,
GOSAT, ...



Research Center Karlsruhe IMK-IFU, Garmisch-Partenkirchen, Ralf Sussmann

Water vapor retrieval from Zugspitze FTIR measurements

Zugspitze/Garmisch instrumentation for water vapor columns/profiles




Research Center Karlsruhe IMK-IFU, Garmisch-Partenkirchen, Ralf Sussmann

Water vapor retrieval from Zugspitze FTIR measurements



Zugspitze operational since 1995
typ. 130-140 measurement
days per year

H₂O columns and profiles



Zugspitze
2964 m

Garmisch operational
since 2004
typ. 130-140 measurement
days per year

“Differential FTIR”
with Zugspitze:
H₂O columns and profiles






Garmisch
734 m



Research Center Karlsruhe IMK-IFU, Garmisch-Partenkirchen

Water vapor retrieval from Zugspitze FTIR measurements

Interference errors and micro window selection: Interference errors (I)




Example CO retrieval:
interfering species
O₃, N₂O, CO₂, H₂O

Research Center Karlsruhe IMK-IFU, Garmisch-Partenkirchen, Ralf Sussmann

Water vapor retrieval from Zugspitze FTIR measurements

Interference errors and micro window selection: Interference errors (II)





Example CO retrieval

taken from:
Sussmann, R. and Borsdorff, T:
Technical note: Interference errors in
infrared remote sounding of the
atmosphere,
Atmos. Chem. Phys., 7, 3537-3557,
2007.


Research Center Karlsruhe IMK-IFU, Garmisch-Partenkirchen, Ralf Sussmann

Water vapor retrieval from Zugspitze FTIR measurements

Zugspitze FTIR water vapor profile retrieval: Interference-free micro windows




HITRAN2000



Research Center Karlsruhe IMK-IFU, Garmisch-Partenkirchen, Ralf Sussmann

Water vapor retrieval from Zugspitze FTIR measurements


Zugspitze FTIR water vapor profile retrieval: A priori information used (I)



Research Center Karlsruhe IMK-IFU, Garmisch-Partenkirchen, Ralf Sussmann

Water vapor retrieval from Zugspitze FTIR measurements


Zugspitze FTIR water vapor profile retrieval: A priori information used (II)



Research Center Karlsruhe IMK-IFU, Garmisch-Partenkirchen, Ralf Sussmann

Water vapor retrieval from Zugspitze FTIR measurements


Zugspitze FTIR water vapor profile retrieval: Averaging kernels



Research Center Karlsruhe IMK-IFU, Garmisch-Partenkirchen, Ralf Sussmann

Water vapor retrieval from Zugspitze FTIR measurements

Zugspitze FTIR water vapor profile retrieval: Smoothing error




⇒ Smoothing error/altitude range does not depend upon absolute column level

Research Center Karlsruhe IMK-IFU, Garmisch-Partenkirchen, Ralf Sussmann

Water vapor retrieval from Zugspitze FTIR measurements

Zugspitze FTIR water vapor profile retrieval: Retrieved profile versus sonde



Research Center Karlsruhe IMK-IFU, Garmisch-Partenkirchen, Ralf Sussmann

Water vapor retrieval from Zugspitze FTIR measurements

Strategy for validation of FTIR total water vapor columns by sondes: Tobin-Sondes

Sonde 1 launched 1h before overpass
Sonde 2 launched 5 min before overpass

AIRS validation campaign
19 Aug 2002 - 17 Nov 2002

Vaisala RS 80-30 G sondes
TOTEX-800-g balloons
2 x Digicora III (Marvin 21, SPS220G)



TOBIN-Inter-/Extrapolation between both soundings:

$$q_{\text{Tobin}}(z, t_{\text{op}}) = q_{\text{sonde}}(z, t_0) + (dq(z)/dt) (t_{\text{op}} - t_0)$$


Tobin, D., W. Feltz, B. Knuteson, H. Revercomb, "ARM T/q Best Estimate Profiles for AIRS validation", 1 March 2000



Research Center Karlsruhe IMK-IFU, Garmisch-Partenkirchen

Water vapor retrieval from Zugspitze FTIR measurements


Strategy for validation of FTIR total water vapor columns by sondes: [Tobin-Sondes](#)



Research Center Karlsruhe IMK-IFU, Garmisch-Partenkirchen, Ralf Sussmann

Water vapor retrieval from Zugspitze FTIR measurements

Strategy for validation of FTIR total water vapor columns by sondes: [Tobin-Sondes](#)




Research Center Karlsruhe IMK-IFU, Garmisch-Partenkirchen, Ralf Sussmann

Water vapor retrieval from Zugspitze FTIR measurements

Optimized FTIR total water vapor column retrieval and its validation

Columns
above
Zugspitze,
2964 m

2-h-mean
values




⇒ FTIR is one of the most accurate techniques to quantify columnar water vapor

Research Center Karlsruhe IMK-IFU, Garmisch-Partenkirchen, Ralf Sussmann

Water vapor retrieval from Zugspitze FTIR measurements

Optimized FTIR total water vapor column retrieval and its validation



Research Center Karlsruhe IMK-IFU, Garmisch-Partenkirchen, Ralf Sussmann

Water vapor retrieval from Zugspitze FTIR measurements

Summary

- found interference-free micro-window set
- profile retrieval
 - OE with climatological a priori covariance from sondes: dofs ≈ 3 , i.e., layers 3 - 4 km, 4 - 7 km, 7 - 15 km
 - smoothing error/altitude range does not depend on absolute water column
 - FTIR profiles reproduce sonde fine structure well, given high quality spectra (snr > 120) and dry days (< 3 mm)
- FTIR is the most accurate technique for total columns (?): validation against sondes (Tobin): $R^2 = 0.98$, slope = 1.01, stdv < 0.3 mm, intercept < 0.01 mm, bias < 0.05 mm

Outlook

- validation against lidar / synergetic combination with lidar
- investigate “differential FTIR“ (Garmisch – Zugspitze column)

IMK-IFU, Research Center Karlsruhe, Garmisch-Partenkirchen, Ralf Sussmann

Water vapor retrieval from Zugspitze FTIR measurements